

DISINFECTION PROCESS

Quick Reference Guide

What You'll Need:

- KlorKleen tablets
- Sterile water
- Gray loop back hose(s)
- Drain pan



STEP 1: DRAIN THE SYSTEM

- Power off and position over drain pan
- Open drain plug
- Unplug one side of the gray loop back hose
- Allow system to fully drain



STEP 2: DISINFECTION CYCLE

- Add 2 KlorKleen tablets to reservoir
- Fill with ~3.8L (1 gallon) sterile water
- Manual Mode: Set temperature to 25°C (77°F)
- Run for a minimum of 20 minutes



STEP 3: RINSE SYSTEM

- Drain system completely
- Refill with ~3.8L sterile water
- Manual Mode: Set temperature to 25°C (77°F)
- Run for a minimum of 5 minutes



STEP 4: FINISH & STORE

- Drain system
- Store Unit Dry

IMPORTANT

Disinfect internal water circuit
and hoses every 14 days

Scan QR code for full
instructions and resources

