

# Altrix Precision Temperature Management System

## Operations Manual

**REF** 8001




EN




# Global symbol glossary

See the Global Symbol Glossary at [3tmedical.com/technical-docs](http://3tmedical.com/technical-docs) for symbol definitions.

## Symbols

|                                                                                     |                                                                                                                                                                                                                                       |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Refer to instruction manual/booklet                                                                                                                                                                                                   |
|    | General mandatory action sign                                                                                                                                                                                                         |
|    | Consult instructions for use                                                                                                                                                                                                          |
|    | General warning                                                                                                                                                                                                                       |
|    | Caution                                                                                                                                                                                                                               |
|    | Warning; electricity                                                                                                                                                                                                                  |
|    | Catalogue number/model                                                                                                                                                                                                                |
|   | Serial number                                                                                                                                                                                                                         |
|  | Unique device identifier                                                                                                                                                                                                              |
|  | Quantity                                                                                                                                                                                                                              |
|  | For US Patents see <a href="http://3tmedical.com/patents">3tmedical.com/patents</a>                                                                                                                                                   |
|  | Manufacturer                                                                                                                                                                                                                          |
|  | Date of manufacture                                                                                                                                                                                                                   |
|  | Mass of equipment                                                                                                                                                                                                                     |
|  | Direct current                                                                                                                                                                                                                        |
|  | Alternating current                                                                                                                                                                                                                   |
|  | Product provides terminal for connection of a potential equalization conductor. The potential equalization conductor provides direct connection between the product and potential equalization busbar of the electrical installation. |
|  | Protective earth ground                                                                                                                                                                                                               |
| IPX1                                                                                | Protection from dripping water from above the device                                                                                                                                                                                  |

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Defibrillation proof type BF (Body Floating) applied part                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|    | <b>CAUTION</b> - Federal law (USA) restricts this device to sale by or on the order of a physician.                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|    | <b>CAUTION</b> - Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.                                                                                                                                                                                                                                                                                                                                                                             |
|    | In accordance with European Directive 2012/19/EU on Waste Electrical and Electronic Equipment (WEEE) as amended, this symbol indicates that the product should be collected separately for recycling. Do not dispose of as unsorted municipal waste. Contact local distributor for disposal information. Ensure infected equipment is decontaminated prior to recycling.                                                                                                                                                                       |
|    | Medical Equipment Classified by Underwriters Laboratories Inc. with Respect to Electric Shock, Fire, Mechanical and Other Specified Hazards Only in Accordance with IEC 60601-1:20 05 (3rd edition), ANSI/AAMI ES60601-1 (2005, 3rd edition), CAN/CSA C22.2 No. 60601-1:20 08, IEC 80601-2-35:2009, CAN/CSA C22.2 NO 80601-2-35:12, ISO 80601-2-56:2009, CAN/CSA C22.2 NO 80601-2-56:12, IEC 60601-1-8:2007, CAN/CSA C22.2 NO 60601-1-8-08, IEC 60601-1-10:2008, CAN/CSA C22.2 NO 60601-1-10-09, IEC 60601-1-6, CAN/CSA-C22.2 No. 60601-1-6:11 |
|    | Liquid level indicator                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|  | Fragile, handle with care                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |
|  | Keep dry                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|  | Do not stack                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|  | This way up                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |

# Table of Contents

|                                                                |    |
|----------------------------------------------------------------|----|
| Warning/Caution/Note Definition .....                          | 3  |
| Summary of safety precautions .....                            | 3  |
| Introduction .....                                             | 5  |
| Product description .....                                      | 5  |
| Indications for use .....                                      | 6  |
| Intended users .....                                           | 6  |
| Clinical benefits .....                                        | 6  |
| Contraindications .....                                        | 6  |
| Expected service life .....                                    | 6  |
| Disposal/recycle .....                                         | 6  |
| Specifications .....                                           | 6  |
| Refrigerant specifications .....                               | 8  |
| Product illustration .....                                     | 9  |
| Product system .....                                           | 10 |
| Product functions .....                                        | 11 |
| Buttons .....                                                  | 11 |
| Visual indicators .....                                        | 12 |
| Flow indicators .....                                          | 13 |
| Graphical user interface icons .....                           | 13 |
| Product alarms .....                                           | 14 |
| Alarm priority and description .....                           | 14 |
| Contact information .....                                      | 16 |
| Serial number location .....                                   | 17 |
| Date of manufacture .....                                      | 17 |
| Setup .....                                                    | 18 |
| Inspecting .....                                               | 18 |
| Selecting a language .....                                     | 18 |
| Testing visual and audible alarms .....                        | 18 |
| Operation .....                                                | 20 |
| Placing the product .....                                      | 20 |
| Applying or releasing the wheel locks .....                    | 20 |
| Selecting and connecting a temperature probe .....             | 20 |
| Connecting the reusable patient temperature output cable ..... | 21 |
| Connecting the insulated hoses .....                           | 22 |
| Disconnecting the insulated hoses .....                        | 22 |
| Connecting and disconnecting thermal transfer devices .....    | 22 |
| Powering on the product .....                                  | 24 |
| Removing and replacing the reservoir .....                     | 24 |
| Filling the reservoir with sterile distilled water .....       | 25 |
| Selecting and setting the primary probe .....                  | 26 |
| Filling a thermal transfer device .....                        | 26 |
| Selecting a therapy mode .....                                 | 27 |
| Starting Automatic therapy mode .....                          | 27 |
| Setting or editing the cooling rates .....                     | 28 |
| Setting or editing the warming rates .....                     | 28 |
| Starting Manual mode .....                                     | 29 |
| Starting Monitor mode .....                                    | 29 |
| Switching modes .....                                          | 30 |
| Pausing and resuming therapy .....                             | 30 |
| Displaying the data storage .....                              | 30 |
| Opening and securing items in the storage compartment .....    | 31 |
| Stopping therapy or powering off the product .....             | 31 |
| Draining the thermal transfer devices .....                    | 31 |
| Draining water from the reservoir .....                        | 32 |
| Draining water from the controller and hoses .....             | 32 |
| Emptying water with a drain pan .....                          | 33 |

|                                                                       |    |
|-----------------------------------------------------------------------|----|
| Storing the power cord and hoses .....                                | 34 |
| Storing the controller .....                                          | 34 |
| Transporting the product.....                                         | 35 |
| Accessories and parts .....                                           | 36 |
| Thermal transfer devices .....                                        | 36 |
| Thermal transfer device kits .....                                    | 36 |
| Patient temperature probes .....                                      | 37 |
| Cables .....                                                          | 37 |
| Hoses .....                                                           | 37 |
| Cleaning tools.....                                                   | 38 |
| Preventive maintenance .....                                          | 39 |
| Cleaning.....                                                         | 40 |
| Cleaning the external surfaces .....                                  | 40 |
| Disinfecting .....                                                    | 41 |
| Disinfecting external surfaces .....                                  | 41 |
| Disinfecting the internal water circuit and hoses every 14 days ..... | 43 |
| Draining the internal water circuit and hoses for disinfection .....  | 43 |
| Disinfecting the internal water circuit and hoses.....                | 45 |
| Rinsing the internal water circuit and hoses.....                     | 47 |
| Troubleshooting .....                                                 | 49 |
| Alarm conditions .....                                                | 51 |
| Check patient probe alarm.....                                        | 52 |
| Patient probe malfunction alarm .....                                 | 52 |
| Patient probe disconnect alarm .....                                  | 52 |
| Patient temperature deviation medium alarm .....                      | 52 |
| Patient temperature deviation low alarm .....                         | 53 |
| Patient temperature output deviation alarm .....                      | 53 |
| Normothermia deviation alarm .....                                    | 53 |
| Water temperature deviation alarm group .....                         | 53 |
| Check water flow alarm .....                                          | 53 |
| No water flow alarm .....                                             | 54 |
| No water alarm .....                                                  | 54 |
| Power backup level alarm.....                                         | 54 |
| Therapy paused time out alarm .....                                   | 55 |
| Remove from use mode .....                                            | 55 |
| EMC Information.....                                                  | 56 |

# Warning/Caution/Note Definition

The words **WARNING**, **CAUTION**, and **NOTE** carry special meanings and should be carefully reviewed.

## **WARNING**

Alerts the reader about a situation which, if not avoided, could result in death or serious injury. It may also describe potential serious adverse reactions and safety hazards.

## **CAUTION**

Alerts the reader of a potentially hazardous situation which, if not avoided, may result in minor or moderate injury to the user or patient or damage to the product or other property. This includes special care necessary for the safe and effective use of the device and the care necessary to avoid damage to a device that may occur as a result of use or misuse.

**Note** - Provides special information to make maintenance easier or important instructions clearer.

## Summary of safety precautions

Always read and strictly follow the warnings and cautions listed on this page. Service only by qualified personnel.

---

### **WARNING**

- Always use 3T Medical accessories. Only IEC 60601-1 equipment shall be hooked to the patient temperature ports. Failure to comply with these instructions may invalidate any or all warranties and may negatively affect the product's EMC performance. This also protects the product from cardiac defibrillation.
  - Avoid reduction in water flow. Do not connect two or more thermal transfer devices in a series on a single port.
  - When you operate the product near ambient temperature limitations of 15.0° C (59.0° F) or 32.0° C (89.6° F), you may experience a reduction in product performance.
  - Always turn or re-position the patient over the duration of therapy, if possible, to reduce the risk of pressure ulcers. Follow your hospital protocol.
  - Always check the integrity of the patient's skin and temperature according to hospital protocol when using the Altrix system.
  - When using the temperature controlled Automatic therapy mode for warming (min, med, or custom), switching to other modes, changing the target patient temperature, or changing the therapy selection may impact the overall benefit of therapy.
- 

### **CAUTION**

- Federal law (USA) restricts this device to sale by or on the order of a physician.
- Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.
- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
- Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- Shock Hazard. Improper handling of the power cord may damage the power cord and cause potential shock hazards. If damage has occurred to the power cord, immediately remove the temperature management system from service to avoid the risk of serious injury or death. Contact the appropriate maintenance personnel.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment like Altrix. Install and place Altrix into service according to the EMC information located in the EMC section of this manual. Portable and mobile RF communications equipment can affect the function of Altrix.
- Shock Hazard. If the internal electrical components are exposed, because the side panel or cover are compromised, remove the product from use.
- Always make sure that the product reaches room temperature before you setup or operate the product.
- Before first use, disinfect the internal water circuit.
- Do not use Altrix located near or stacked with other medical equipment. If it is necessary to locate Altrix near other medical equipment, make sure it operates as intended.

- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be used no closer than 12 inches (30 cm) to any part of Altrix system, including cables specified by the manufacturer.
  - Always apply the wheel locks to prevent unintended movement.
  - Do not use high frequency surgical instruments or endocardial catheters while the Altrix system is in use. This is to avoid the risk of electrical shock, burns, or electromagnetic interference.
  - Avoid the use of materials of good thermal conductivity, such as water, gel, or similar substances, with the Altrix system powered off. This can decrease the temperature of the body of a patient.
  - Do not apply thermal transfer devices to patients with ischemic limbs. This may result in harm to the patient.
  - Do not use this product if the patient has a transdermal medication (patch) as this can result in increased drug delivery.
  - Do not use three or more adult Mul-T-Blanket products at the same time to avoid the risk of water overflow when you power off the controller.
  - Always prefill the thermal transfer devices with sterile, distilled water before you apply it to the patients. Prefilling reduces the risk of pressure ulcers.
  - Always clamp the hoses when you disconnect the thermal transfer devices.
  - Electric shock. This equipment must only be connected to a supply mains with protective earth.
  - Always plug this product directly into a properly grounded hospital-grade or medical-grade wall outlet to achieve grounding reliability.
  - Do not use high frequency surgical instruments or endocardial catheters while the Altrix system is in use. This is to avoid the risk of electrical shock, burns, or electromagnetic interference.
  - Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal or mask type.
  - Do not place cables, hoses, or power cord in walkways to avoid the risk of trip hazards.
  - Do not place your fingers in between the reservoir and the sides of the controller, to avoid the risk of pinching your fingers.
  - Always fill the reservoir with room temperature sterile distilled water to reduce the risk of burn.
  - Do not overfill the reservoir to avoid the risk of water spillage and fall.
  - Always make sure that there are no water leaks before starting a defibrillation.
  - Always monitor the patient for shivering, temperature, signs of intolerance, and skin condition when using this product.
  - Do not hang items on the controller handle to avoid the risk of tipping the product.
  - Always store the power cord, cables, and hoses before you transport the product to reduce the risk of a trip hazard.
  - Do not store the product with water in the device.
  - Always store the product within the specified environmental condition values.
  - Always use extra care when you transport the product long distances and on inclines greater than five degrees. Ask for help, if necessary, to avoid the risk of tipping.
  - Always use the handle to move the product. Do not attempt to move the product by pulling on cables, hoses, or by any other means.
  - Avoid ramps that are steeper than ten degrees to avoid tipping the product.
  - Do not power wash this product.
  - Do not use quaternaries that contain glycol ethers as they may damage the reusable accessories.
  - Do not disinfect the internal water system with a thermal transfer device attached as this may cause a leak.
  - Do not use bleach or any other cleaning or disinfectant agents for internal circuits. This could result in damage to the product. Only use approved disinfectant tablets.
  - Always drain the product before disinfecting the internal water circuit. Failure to drain the product may reduce the effectiveness of the disinfection process.
  - Always remove the product from use before servicing any components. Contact qualified service personnel for service.
  - This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.
- 

**Note** - Disinfection of the Altrix internal water system was validated using *M. mucogenicum*.

# Introduction

This manual assists you with the operation or maintenance of your 3T Medical product. Read this manual before operating or maintaining this product. Set methods and procedures to educate and train your staff on the safe operation or maintenance of this product.

---

## CAUTION

- Improper usage of the product can cause injury to the patient or operator. Operate the product only as described in this manual.
  - Do not modify the product or any components of the product. Modifying the product can cause unpredictable operation resulting in injury to patient or operator. Modifying the product also voids its warranty.
- 

## Note

- This manual is a permanent part of the product and should remain with the product even if the product is sold.
- 3T Medical continually seeks advancements in product design and quality. This manual contains the most current product information available at the time of printing. There may be minor discrepancies between your product and this manual. If you have any questions, contact 3T Medical Customer Service or Technical Support at (734) 525-9580.

## Product description

The 3T Medical Model 8001 Altrix Precision Temperature Management System can supply water to an individual or multiple thermal transfer devices simultaneously with each of these circuits monitored separately. Three operating modes are available to ease patient care: Automatic, Manual, and Monitor. The controller uses the patient temperature probe to provide closed loop feedback for automatic patient temperature management and monitoring. The controller alarms visual and audible indications for when safety parameters are exceeded or it detects system function or performance irregularities. The Altrix system is able to provide a patient temperature output reference signal to be connected to a non-specific third party device or system.

The controller regulates water temperatures between 4.0° C (39.2° F) and 40.0° C (104.0° F) and circulates the heated or cooled water via hose sets through the thermal transfer devices. A graphical display provides the user an interface for selecting desired water or patient temperature settings, operating modes, help menus, and other key parameters. Visual indicators are displayed to inform the user of system status or when the user must confirm a setting selection. The system's water temperature and flow outputs can be monitored with 400 series compatible devices to optimize system operation.

The Altrix system includes the following components:

- controller
- reusable hose sets
- thermal transfer devices (blankets, vests, and leg wraps)
- patient temperature probes
- reusable adapter cables
- reusable patient temperature output cable

## Note

- The blankets, vests, leg wraps, and patient temperature probes are type BF applied parts.
- Sterile distilled water = Sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns.

## Indications for use

The Altrix system is intended for circulating temperature controlled warm or cold water via patient contact thermal transfer devices for the application of regulating human body temperature in situations where a physician or clinician with prescription privileges determines that temperature therapy is necessary or desirable.

Indications for use for the Altrix system include:

- Maintain pre-set body temperature as determined by the physician
- Maintain normal body temperature during surgical procedures
- For use in all clinical settings including coronary care units, operating, recovery and emergency departments, burn units, and medical/surgical units
- Adult and pediatric patients
- Monitoring and controlling patient temperature
- Temperature reduction in patients where clinically indicated, e.g., in hyperthermic patients

## Intended users

- Physicians
- Advanced Practice Registered Nurses
- Nurses

## Clinical benefits

Regulating patient temperature

## Contraindications

For core body temperature regulation:

- Raynaud's Phenomenon (primary or secondary)
- Application to lower extremities distal to aortic cross-clamping

## Expected service life

The Altrix controller has a five year expected service life under normal use conditions and with appropriate periodic maintenance. See the maintenance manual for preventive maintenance and service information.

## Disposal/recycle

Always follow the current local recommendations and/or regulations governing environmental protection and the risks associated with recycling or disposing of the equipment at the end of its useful life.

## Specifications









|                                                                        |                     |
|------------------------------------------------------------------------|---------------------|
| <b>Model</b>                                                           | <b>8001-000-001</b> |
| Electrical Requirements - AC Voltage Input Current and Voltage Ratings | 120VAC, 60Hz 12A    |

| Physical dimensions               |                                                                                                                                       |                  |
|-----------------------------------|---------------------------------------------------------------------------------------------------------------------------------------|------------------|
| Height                            | 42.5 in.                                                                                                                              | 107.9 cm         |
| Width                             | 15.0 in.                                                                                                                              | 38.1 cm          |
| Depth                             | 23.0 in.                                                                                                                              | 58.4 cm          |
| Empty weight                      | 150.0 lb                                                                                                                              | 68.0 kg          |
| Filled weight                     | 160.5 lb                                                                                                                              | 72.8 kg          |
| Reservoir capacity                | 1.3 gal                                                                                                                               | 5.0 L            |
| Water temperature                 |                                                                                                                                       |                  |
| Control setting range             | 39.2° - 104.0° F                                                                                                                      | 4.0° - 40.0° C   |
| Control accuracy                  | ±0.3° C (4.0° - 40.0° C)                                                                                                              |                  |
| Display measurement accuracy      | ±0.2° C (4.0° - 40.0° C)                                                                                                              |                  |
| Display/resolution setting        | 0.1° C                                                                                                                                |                  |
| Default setting                   | 104.0° F                                                                                                                              | 40.0° C          |
| Patient temperature               |                                                                                                                                       |                  |
| Control setting range             | 89.6° - 100.4° F                                                                                                                      | 32.0° - 38.0° C  |
| Control accuracy                  | ±0.1° C (32.0° - 38.0°C)                                                                                                              |                  |
| Measurement accuracy              | ±0.3° C (25.0° - 45.0° C)                                                                                                             |                  |
|                                   | ±0.4° C (0° - 24.9° C, 45.1° - 50.0° C)                                                                                               |                  |
| Display/resolution setting        | 0.1° C                                                                                                                                |                  |
| Display range                     | 32.0° - 122.0° F                                                                                                                      | 0.0° - 50.0° C   |
| Default setting                   | 98.6° F                                                                                                                               | 37.0° C          |
| Controller                        |                                                                                                                                       |                  |
| Heater capacity, maximum          | 500 watts                                                                                                                             |                  |
| Circulating fluid                 | Sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product |                  |
| Battery                           | 9 V Lithium                                                                                                                           |                  |
| Alarm tone range                  | 75 - 85 dBA per standard IEC 60601-1-8                                                                                                |                  |
| Water flow rate in each hose port | Typical 1.2 lpm                                                                                                                       |                  |
| Refrigerant type                  | Units produced before May 2025: R-134a                                                                                                |                  |
|                                   | Units produced after May 2025 or thermal units replaced after May 2025: R-513a                                                        |                  |
| Power cord length                 | 14.0 - 15.0 feet                                                                                                                      | 4.2 - 4.5 meters |
| Clinical thermometer              | Direct mode                                                                                                                           |                  |
| Equipment Class                   | Class I                                                                                                                               |                  |
|                                   | Rate for continuous operation                                                                                                         |                  |

**Note** - The controller takes approximately 9 minutes to heat from 23.0±2.0° C (73.4° F) to 37.0° C (98.6° F) when not connected to a patient. Time will vary when connected to a patient.

3T Medical reserves the right to change specifications without notice.

For more information about thermal transfer devices, cables, or probes, see the manufacturer's instructions for use.

| Environmental conditions           | Operation                                                                         | Storage                                                                            | Transportation                                                                      |
|------------------------------------|-----------------------------------------------------------------------------------|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| Ambient temperature                |  |  |  |
| Relative humidity (non-condensing) |  |  |  |
| Atmospheric pressure               |  |  | Not applicable                                                                      |

## Refrigerant specifications

- Altrix contains a fluorinated greenhouse gas, and its functioning relies upon this gas (refrigerant).
- The industry designation for the fluorinated greenhouse gas used is R-513a.
- The unit is filled with 6.0 - 9.0 ozf (force) refrigerant R-513a into the system. This equates to approximately 0.45 lb (0.20 kg) of refrigerant.
- The global warming potential (GWP) of R-513a is 573.




Figure 1 – Refrigerant label location

**Note** - For medical applications only.

# Product illustration




Figure 2 – Controller, patient front

|   |                                 |
|---|---------------------------------|
| A | Storage compartment             |
| B | Removable water reservoir       |
| C | Front wheel                     |
| D | Patient probe port              |
| E | Patient temperature output port |
| F | Hose connection ports/collars   |
|   |                                 |




Figure 3 – Controller, patient back

|   |                                      |
|---|--------------------------------------|
| G | Graphical user interface display     |
| H | Handle                               |
| I | Power cord                           |
| J | Hose and power cord management strap |
| K | Swivel caster                        |
| L | Wheel lock                           |
| M | Ground chain                         |

Product system



Figure 4 – Altrix system - controller with thermal transfer devices





Figure 5 – Closed loop system

|   |                                      |    |                        |
|---|--------------------------------------|----|------------------------|
| 1 | Human machine interface (HMI) system | 6  | Controls               |
| 2 | Physical boundary                    | 7  | Flow system            |
| 3 | Patient system                       | 8  | Energy transfer system |
| 4 | Patient temperature port             | 9  | Water flow             |
| 5 | Fluid delivery system                | 10 | Signals                |

## Product functions

The graphical user interface shown is for reference only. The image shows where you will see the icons and buttons illuminate when they are active. At no time will you see all of these icons at the same time.



## Buttons

The buttons are located on the outside of the graphical user interface. They are visible when available.

| Icon                                                                                | Name                   | Function                                                                                                                                                                                                                 |
|-------------------------------------------------------------------------------------|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|   | Stand-by               | Press and hold the button for two seconds to stop therapy or power off                                                                                                                                                   |
|  | Therapy paused         | Press and hold the button for two seconds to pause or resume therapy                                                                                                                                                     |
|  | View temperature       | Select temperature degree in Celsius or Fahrenheit                                                                                                                                                                       |
|  | Lock/unlock screen     | Press and hold the button for two seconds to lock or unlock the graphical user interface                                                                                                                                 |
|  | Audio paused           | Pause or resume the audible indicator when an alarm is active. Silences each alarm for 5, 10 or 30 minutes depending on the alarm condition. This button breathes <sup>1</sup> to indicate that it is in a paused state. |
|  | Automatic therapy mode | Cools or warms the patient to a selected patient target temperature                                                                                                                                                      |
|  | Manual therapy mode    | Cools or warms the water to a selected water target temperature                                                                                                                                                          |
|  | Monitor only mode      | Displays the current patient temperature (no therapy)                                                                                                                                                                    |
|  | Increase               | Increases the water or patient temperature by 0.1° for cooling or warming temperature<br><b>Note</b> - Press and hold the increase button to move the temperature up faster.                                             |

| Icon                                                                              | Name                                       | Function                                                                                                                                                                       |
|-----------------------------------------------------------------------------------|--------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Decrease                                   | Decreases the water or patient temperature by 0.1° for cooling or warming temperature<br><b>Note</b> - Press and hold the decrease button to move the temperature down faster. |
|  | Back                                       | Returns to the previous screen or cancel an operation                                                                                                                          |
|  | Edit settings, Exit, or Cancel             | Edit current settings, exit, or cancel                                                                                                                                         |
|  | Confirm selection                          | Accepts the selected settings                                                                                                                                                  |
|  | Next or More                               | Changes to the next screen, option, or setting                                                                                                                                 |
|  | Page indicators (may also appear vertical) | Indicates that there is more than one page associated with the screen topic for the page that is currently displayed                                                           |
|  | Settings                                   | Displays the summary of the current, visual/audible, language, or primary probe settings                                                                                       |
|  | Graph                                      | Graphical display of the selected items such as patient temperature, target temperature, water temperature, and power level                                                    |
|  | Help                                       | Displays contextual help screens for therapies, navigation, buttons, and alarm screens. This button breathes to allow the user to view the alarm screen.                       |

#### Note

- If not specified above, make sure that you press and release the buttons or icons for your selection to register with the system.
- The Light sensor (non selectable)  dims or brightens the LCD based on the amount of light in the room.
- <sup>1</sup>Breathe: The brightness of the button or icon will go to a low light and then increase to a bright light. This cycle repeats.

#### Visual indicators

When the visual indicators are solid green, this indicates that the function is stabilized. The visual indicators breathe to indicate that the controller is at the intermediate target.



| Icon, green                                                                         | Name                                                                   |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|  | Water temperature on target, solid green when active, does not breathe |
|  | Patient temperature on target                                          |
|  | Patient probe A port, stabilized                                       |
|  | Patient probe B port, stabilized                                       |




| Icon, green                                                                       | Name                                                              |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------|
|  | External device, patient probe A                                  |
|  | External device, patient probe B                                  |
|  | Stand-by                                                          |
|  | Water flow detected, ports 1, 2, or 3 are solid green when active |
|                                                                                   | If the water flow is less than 0.8 lpm, ports 1, 2, or 3 breathe  |

## Flow indicators

| Indication               | Manual Mode                                                                                                                                                                                                                                                                                      | Automatic Mode                                                                                                                                                                                                                                              |
|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Low flow <0.8 lpm        | Indicated by a breathing green indicator for 10 minutes, maximum<br><br>Low flow that changes to No flow will remain breathing green for 15 seconds, maximum                                                                                                                                     | Indicated by a breathing green indicator (no alarms)<br><br>Low flow that changes to No flow will remain breathing green for 15 seconds, maximum (no alarms)                                                                                                |
| Good flow $\geq$ 0.8 lpm | Solid green indicator when flow is maintained for 30 seconds.<br><br>Good flow that is reduced to Low flow will change the indicator to breathing green for 10 minutes, maximum<br><br>Good flow that is reduced to No flow will change the indicator to breathing green for 15 seconds, maximum | Solid green indicator when flow is maintained for 30 seconds.<br><br>Good flow that is reduced to Low flow (no alarm)<br><br>Good flow that is reduced to No flow and will change the indicator to breathing green after 15 seconds the indicator turns off |

## Graphical user interface icons

| Icon                                                                                | Name                     |
|-------------------------------------------------------------------------------------|--------------------------|
|  | Cooling therapy          |
|  | Warming therapy          |
|  | Current therapy duration |
|  | Total duration           |
|  | Visual and audible tests |

| Icon                                                                              | Name                                                                                                                                                              |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | Target patient or water temperature                                                                                                                               |
|  | Medium: patient temperature increases at a rate of 4.0° C in 12 hours (0.33° C/ hour).                                                                            |
|  | Maximum: water temperature approaches water target as fast as possible                                                                                            |
|  | Minimum: patient temperature increases at a rate of 4.0° C in 24 hours (0.17° C/hour).                                                                            |
|  | Set Custom: patient temperature increases at a customized temperature and time period the operator selects. The temperature increases 0.05° C/hour to 0.5° C/hour |
|  | Maximum: water temperature approaches water target as fast as possible                                                                                            |
|  | Medium: water is cooled to target, with a max of 15.0° C difference between the patient and the water temperature                                                 |
|  | Minimum: water is cooled to target, with a max of 10.0° C difference between the patient and the water temperature                                                |

## Product alarms

Audible alarms work in conjunction with the display.

### Alarm priority and description

| Priority alarm | Audible reminder                                                                                                                                                                                                                                                                                                                 | Icon flashes                                                                                                                                                     |
|----------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Medium         | <p>Repetitive burst of 3 beeps upon first occurrence and every 25 seconds until the condition is resolved.</p> <p>If you pause the audio, when the timer expires, the audio indication will resume.</p>                                                                                                                          | When a medium priority alarms, the icon will breathe yellow to indicate there is an alarm. The yellow icon will continue to breathe until the alarm is resolved. |
| Low            | Single burst of 2 beeps upon first occurrence.                                                                                                                                                                                                                                                                                   | The icon for the low priority alarm is solid yellow. The solid yellow icon will remain until the alarm condition is resolved.                                    |
| Audio pause    | <p>The audio pause button breathes as a reminder and resets all of the active alarms.</p> <p>The audible alarm will resume as soon as a new alarm is raised or the last Medium audio paused timer has expired.</p> <p>You can tap the audio pause button at any time during the alarm condition, as many times as necessary.</p> | Pausing the alarm will not stop the yellow icon from breathing or being solid yellow.                                                                            |

| Icon, yellow                                                                                                                                                           | Name                                  | Alarm priority and delay | Message                                                                                                                                           | Therapy interrupted | Check                                                                                                                                                |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| <br> | Water temperature deviation           | Medium, 5 minute delay   | Water temperature is $\pm 0.8^{\circ}\text{C}$ ( $1.4^{\circ}\text{F}$ ) outside of target temperature                                            | No                  | Temporary condition upon startup, addition of thermal transfer device, or addition of water                                                          |
|                                                                                                                                                                        | No water                              | Medium, 60 second delay  | No water                                                                                                                                          | Yes                 | Check for leaks<br><br>Add a minimum of 2 liters of water                                                                                            |
|                                                                                                                                                                        | No water flow                         | Medium, 20 second delay  | No flow detected                                                                                                                                  | Yes                 | Check for leaks and obstructions at connections, hoses, and thermal transfer devices                                                                 |
|                                                                                                                                                                        | Check water flow                      | Medium, 10 minute delay  | Reduced flow detected                                                                                                                             | No                  | Tap Confirm, if the water port was removed intentionally<br><br>Check for leaks and obstructions at connections, hoses, and thermal transfer devices |
|                                                                                      | Check patient probe (A or B)          | Medium                   | Abnormal change in patient temperature                                                                                                            | Yes                 | Check probe condition, location, and connections                                                                                                     |
|                                                                                                                                                                        | Probe or adapter malfunction (A or B) | Medium, 30 second delay  | No temperature signal detected.                                                                                                                   | Yes                 | Check probe or adapter cable condition, location, and connections                                                                                    |
|                                                                                                                                                                        | Adapter cable disconnected (A or B)   | Medium, 30 second delay  | Adapter cable not detected                                                                                                                        | Yes                 | Reinsert the adapter cable. If damaged, replace the adapter cable                                                                                    |
|                                                                                                                                                                        | Patient temperature deviation         | Medium                   | Patient temperature is more than $\pm 1.0^{\circ}\text{C}$ from target temperature<br><br>(Only will appear after the initial target is reached.) | No                  | Check patient condition, placement of thermal transfer devices, and all connections                                                                  |
|                                                                                                                                                                        | Patient temperature deviation         | Low                      | Patient temperature is more than $\pm 0.5^{\circ}\text{C}$ from target temperature<br><br>(Only will appear after the initial target is reached.) | No                  | Check patient condition, placement of thermal transfer devices, and all connections                                                                  |

| Icon, yellow                                                                        | Name                                | Alarm priority and delay | Message                                                                                     | Therapy interrupted | Check                                                                                                     |
|-------------------------------------------------------------------------------------|-------------------------------------|--------------------------|---------------------------------------------------------------------------------------------|---------------------|-----------------------------------------------------------------------------------------------------------|
|                                                                                     | Normothermia deviation              | Low                      | Patient temperature is outside of 36.0° C (96.8° F) to 38.0° C (100.4° F)                   | No                  | Check patient condition, placement of thermal transfer devices, and all connections                       |
|    | Therapy pause                       | Medium                   | Therapy is currently paused                                                                 | Yes                 | To resume, tap and hold the play/pause for 2 seconds                                                      |
|                                                                                     | Battery low                         | Low                      | Battery is low                                                                              | No                  | Maintenance is recommended. If battery is not replaced, the product may not function on the next startup. |
|                                                                                     | Patient temperature output (A or B) | Low                      | Patient temperature output is inaccurate on the external device, or outside supported range | No                  | Check output adapter cable connection. Tap Confirm to reactivate the output port.                         |
|    | Remove from use (RFU)               | Medium                   | The system has powered off due to a malfunction                                             | Yes                 | Remove the product from use immediately. Notify the appropriate personnel.                                |
|  | Power loss                          | Medium                   | Not applicable                                                                              | Yes                 | Check power cord connection                                                                               |

#### Note

- If any of the alarm conditions persist, call maintenance.
- If page indicators appear on the alarm screen, there are multiple active alarms. The highest alarm displays. Tap Next or Back to view the active alarms.

## Contact information

Contact 3T Medical Customer Service or Technical Support at: (734) 525-9580.


3T Medical Systems, Inc.  
236 Fernwood Avenue  
Edison, NJ 08837  
USA

**Note** - The user and/or the patient should report any serious product-related incident to both the manufacturer and the Competent authority of the European Member State where the user and/or patient is established.

To view your operations or maintenance manual online, see [3tmedical.com/technical-docs](http://3tmedical.com/technical-docs).

Have the serial number (A) of your 3T Medical product available when calling 3T Medical Customer Service or Technical Support. Include the serial number in all written communication.

## Serial number location



## Date of manufacture

The year of manufacture is the first four digits of the serial number.

# Setup

Unpack the cartons and check all items. Make sure that the product is free from visual damage before putting into service.

## CAUTION

- Shock Hazard. Improper handling of the power cord may damage the power cord and cause potential shock hazards. If damage has occurred to the power cord, immediately remove the temperature management system from service to avoid the risk of serious injury or death. Contact the appropriate maintenance personnel.
- Take special precautions regarding electromagnetic compatibility (EMC) when using medical electrical equipment like Altrix. Install and place Altrix into service according to the EMC information located in the EMC section of this manual. Portable and mobile RF communications equipment can affect the function of Altrix.
- Shock Hazard. If the internal electrical components are exposed, because the side panel or cover are compromised, remove the product from use.
- Always make sure that the product reaches room temperature before you setup or operate the product.
- Before first use, disinfect the internal water circuit.

## Inspecting

Before you place the product into service, make sure that the controller works.

1. Visually inspect the product for any signs of shipping damage.
2. Plug the power cord into a wall outlet. Make sure that the power indicator illuminates on the operator control panel.
3. Before first use, see *Disinfecting the internal water circuit and hoses every 14 days* (page 43).

## Selecting a language

The Altrix controller has several language choices. English is the default language.



To select a language when in standby mode:

1. Tap the Settings button, to show the Select Language screen.
  - a. Tap Next, if you are in therapy mode.
2. Tap More to view other languages.
3. Select a language. Tap on the Increase or Decrease buttons or tap the language, to highlight a language of your choice.
4. Tap Confirm.

**Note** - If you do not touch the screen for three minutes, the LCD will return to the previous menu.

## Testing visual and audible alarms

Before you place the product into service, make sure that the visual and audible alarms function.

1. Tap the Settings button.
2. Tap the Back button.
3. Tap the Visual/Audible icon.



4. Tap Confirm.

**Note**

- The system runs through visual tests of the Green Indicators, Yellow Indicators, White Indicators, and Fluid Controller Light tests, and audible alarms.
  - The test will continue until you stop it.
5. To stop the Visual/Audible test, press the Back button.
  6. To exit settings, press the Exit button.

# Operation

## Placing the product

When placing the product, do not block access to the hospital-grade plug or medical-grade wall outlet.

### CAUTION

- Do not use Altrix located near or stacked with other medical equipment. If it is necessary to locate Altrix near other medical equipment, make sure it operates as intended.
- Portable RF communications equipment, including peripherals such as antenna cables and external antennas, should be used no closer than 12 inches (30 cm) to any part of Altrix system, including cables specified by the manufacturer.

Place the Altrix controller outside of the patient environment by 1.5 m (4.9 ft) (Figure 6).




Figure 6 – Product placement

## Applying or releasing the wheel locks

The wheel locks are to help keep the product in place. The wheel lock prevents the rear caster wheels from rotating but does not prevent the product from sliding on the floor surface.

**CAUTION** - Always apply the wheel locks to prevent unintended movement.

To apply the wheel locks, push down (A) (Figure 7) with your foot.

To release the wheel locks, pull up (A) (Figure 7) with your foot.




Figure 7 – Wheel locks

## Selecting and connecting a temperature probe


**WARNING** - Always use 3T Medical accessories. Only IEC 60601-1 equipment shall be hooked to the patient temperature ports. Failure to comply with these instructions may invalidate any or all warranties and may negatively affect the product's EMC performance. This also protects the product from cardiac defibrillation.

**CAUTION** - Do not use high frequency surgical instruments or endocardial catheters while the Altrix system is in use. This is to avoid the risk of electrical shock, burns, or electromagnetic interference.

Use only 3T Medical temperature probes. See *Patient temperature probes* (page 37).

To connect the temperature probe:

1. Inspect the temperature probe and reusable adapter cable for wear, breakage, or fraying. Replace, if necessary.
2. Align the red dot on the **Reusable Adapter Cable** (B) to the controller (A) with the red dot on the patient probe port A or port B.



**Figure 8 – Port selected**

3. Connect the plug (C) to the patient temperature probe.
4. Apply the temperature probe to the patient. Follow your hospital protocol and the manufacturer's directions for the selected temperature probe use.
5. Tap Confirm, if applicable.

**Note** - Temperature readings may vary between temperature measurement sites.

## Connecting the reusable patient temperature output cable

This feature allows the operator to view the temperature on the Altrix system and on an external device. Always connect the reusable patient temperature output cable to a 400 series compatible external device for temperature accuracy.

---

**WARNING** - Always use 3T Medical accessories. Only IEC 60601-1 equipment shall be hooked to the patient temperature ports. Failure to comply with these instructions may invalidate any or all warranties and may negatively affect the product's EMC performance. This also protects the product from cardiac defibrillation.

---

To connect the reusable patient temperature output cable:

1. Insert the reusable patient temperature output cable into the patient temperature port (Figure 9).



**Figure 9 – Patient temperature output port**

2. Connect the other end of the reusable patient temperature output cable to the external device.

### **Note**

- When Altrix is powered, the calibration of the patient temperature output is completed.
  - If you need to calibrate the patient temperature output cable, power cycle the product by removing the plug from the wall.
  - For the reusable patient temperature output cable to work properly, make sure that you insert a patient temperature probe into port A or port B.
3. Tap Confirm.

## Connecting the insulated hoses

To connect the insulated hoses:

1. To connect, push back on the retaining collar of the port on the controller (Figure 10).



Figure 10 – Pull back on the retaining collar

2. Push the hose into an upper or lower port (Figure 11) and release the collar until the retaining collar clicks into place (Figure 12).

**Note** - Connect a set of ports for proper water flow.




Figure 11 – Connect the hoses



Figure 12 – Hoses connected

## Disconnecting the insulated hoses

To disconnect the insulated hoses:

1. To disconnect, push back on the retaining collar of the port on the controller.
2. Pull the hose to disconnect.

## Connecting and disconnecting thermal transfer devices

Read the operations manual for the individual thermal transfer devices for warnings, cautions, and safe operating instructions before use.

---

### WARNING

- Always use 3T Medical accessories. Only IEC 60601-1 equipment shall be hooked to the patient temperature ports. Failure to comply with these instructions may invalidate any or all warranties and may negatively affect the product's EMC performance. This also protects the product from cardiac defibrillation.
  - Avoid reduction in water flow. Do not connect two or more thermal transfer devices in a series on a single port.
-

---

## CAUTION

- Avoid the use of materials of good thermal conductivity, such as water, gel, or similar substances, with the Altrix system powered off. This can decrease the temperature of the body of a patient.
  - Do not apply thermal transfer devices to patients with ischemic limbs. This may result in harm to the patient.
  - Do not use this product if the patient has a transdermal medication (patch) as this can result in increased drug delivery.
  - Do not use three or more adult Mul-T-Blanket products at the same time to avoid the risk of water overflow when you power off the controller.
  - Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.
  - Always prefill the thermal transfer devices with sterile, distilled water before you apply it to the patients. Prefilling reduces the risk of pressure ulcers.
  - Always clamp the hoses when you disconnect the thermal transfer devices.
- 

To connect or disconnect the Klik-Tite connectors (Figure 13) to the insulated hoses.




Figure 13 – Klik-Tite

To connect or disconnect the Colder® style (Figure 14) to the insulated hoses.




Figure 14 – Colder style connectors

To close or open hose clamps (Figure 15).

Always clamp the hoses before you disconnect. See *Draining the thermal transfer devices* (page 31).




Figure 15 – Hose clamps

To connect the Altrix Temperature Management Wrap to the Altrix Temperature Management Hose, see Figure 16.

Press the button (A) on the hose to disconnect from the wrap.




Figure 16 – Altrix Temperature Management Wrap connectors

The Altrix Temperature Management Wraps have an internal valve and do not need a clamp to stop the water flow.

**Note** - The term “thermal transfer devices” is used throughout this manual and is interchangeable with blankets and wraps, unless indicated otherwise.

## Powering on the product

The operator should stand in front of the controller within arm's reach. This allows the operator to see and respond to the display notifications.

---

### WARNING

- Avoid reduction in water flow. Do not connect two or more thermal transfer devices in a series on a single port.
  - When you operate the product near ambient temperature limitations of 15.0° C (59.0° F) or 32.0° C (89.6° F), you may experience a reduction in product performance.
- 

### CAUTION

- Shock Hazard. Improper handling of the power cord may damage the power cord and cause potential shock hazards. If damage has occurred to the power cord, immediately remove the temperature management system from service to avoid the risk of serious injury or death. Contact the appropriate maintenance personnel.
  - Electric shock. This equipment must only be connected to a supply mains with protective earth.
  - Always plug this product directly into a properly grounded hospital-grade or medical-grade wall outlet to achieve grounding reliability.
  - Do not use high frequency surgical instruments or endocardial catheters while the Altrix system is in use. This is to avoid the risk of electrical shock, burns, or electromagnetic interference.
  - Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal or mask type.
  - Do not place cables, hoses, or power cord in walkways to avoid the risk of trip hazards.
  - Do not use three or more adult Mul-T-Blanket products at the same time to avoid the risk of water overflow when you power off the controller.
- 

To start the product:

1. Plug the power cord into a wall outlet.



2. Tap the Stand-by button to start the product.
  3. If you are going into Automatic mode or Monitor mode, see *Selecting and setting the primary probe* (page 26). If manual mode, go to the next step.
  4. See *Removing and replacing the reservoir* (page 24).
  5. See *Filling the reservoir with sterile distilled water* (page 25).
  6. Connect up to three thermal transfer devices (with the exception of adult Mul-T-Blankets) to dedicated adapter hoses and ports.
  7. Open the clamps on the connector hose and the thermal transfer devices to provide proper water flow.
  8. See *Filling a thermal transfer device* (page 26).
  9. See *Selecting a therapy mode* (page 27).
  10. Make sure that the desired port configuration is maintained and that water is flowing through the thermal transfer devices.
- 

**WARNING** - Always turn or re-position the patient over the duration of therapy, if possible, to reduce the risk of pressure ulcers. Follow your hospital protocol.

---

## Removing and replacing the reservoir


The removable reservoir enables you to fill or drain the reservoir away from the controller without interrupting therapy. Install the reservoir before you start therapy.

---

**CAUTION** - Do not place your fingers in between the reservoir and the sides of the controller, to avoid the risk of pinching your fingers.

---

To remove the reservoir, pull forward at an angle, and lift out the reservoir (Figure 17).



**Figure 17 – Removable reservoir**

1. To replace the reservoir, align the base of the reservoir over the drain (C).
2. Align the notch on the back of the reservoir (A) with the hook on the controller (B) (Figure 17).
3. Push the reservoir back into place. Make sure that the reservoir is secure to avoid water leakage.

## Filling the reservoir with sterile distilled water

The removable reservoir is translucent for you to see the water levels.

---

### CAUTION

- Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.
  - Always fill the reservoir with room temperature sterile distilled water to reduce the risk of burn.
  - Do not overfill the reservoir to avoid the risk of water spillage and fall.
- 

To fill the removable reservoir with sterile distilled water:

1. See *Removing and replacing the reservoir* (page 24).
2. Fill the reservoir with five liters of sterile distilled water. Do not fill past the top fill line to avoid water overflow (Figure 18).




Figure 18 – Reservoir fill lines

## Selecting and setting the primary probe

A patient probe displays when present, stable, and confirmed. The choice of Probe A or Probe B is highlighted when you insert the cable into port A and port B. If you only insert one cable, the active port displays.

1. Tap the Settings button.
  - a. In the standby mode, tap the Back button to display the edit settings screen.
  - b. In an active therapy mode, tap the Next button.
2. Tap Select Probe to display the Select Primary Probe (Probe A or Probe B) screen. If both probes are present, the default is probe A.
3. Tap A or B, if applicable.
4. Tap Confirm.

### Note

- The message “Probe stabilization in progress... Please wait” is displayed.
- When you initially select a probe (A or B), detected is checked. When stabilized, the Ready check is displayed.
- If the probe is not stabilized within three minutes, the message “Probe stabilization error” will appear. Tap the Help button for more detail.
- You can select help at any time to display help with the current screen or icon descriptions.

## Filling a thermal transfer device

Read the operations manual for the individual thermal transfer devices for warnings, cautions, and safe operating instructions before use.

**Note** - These instructions are to pre-fill the thermal transfer devices only, not therapy. See *Switching modes* (page 30).

To fill a thermal transfer device:

1. Connect a thermal transfer device following: *Connecting and disconnecting thermal transfer devices* (page 22).
2. Lay the thermal transfer device on a flat surface. Make sure the thermal transfer device is flat for water flow.
3. Open all of the clamps on the connector hose and thermal transfer device.
4. Make sure that the controller is powered.
5. Tap the Stand-by button.
6. Tap the Manual mode button.
7. Tap Confirm.

8. Select a water temperature that aligns with your target patient temperature.

**Note** - Allow the water to flow from the controller into the thermal transfer device until full.

9. Tap Confirm.

## Selecting a therapy mode

You can select from one of three therapy modes and tap Confirm:

- Automatic therapy
- Manual therapy
- Monitor non-therapy

For mode descriptions, tap the Help button.

---

### WARNING

- Always check the integrity of the patient's skin and temperature according to hospital protocol when using the Altrix system.
  - Always use 3T Medical accessories. Only IEC 60601-1 equipment shall be hooked to the patient temperature ports. Failure to comply with these instructions may invalidate any or all warranties and may negatively affect the product's EMC performance. This also protects the product from cardiac defibrillation.
  - When using the temperature controlled Automatic therapy mode for warming (min, med, or custom), switching to other modes, changing the target patient temperature, or changing the therapy selection may impact the overall benefit of therapy.
- 

### CAUTION

- Explosion risk. This product is not suitable for use in the presence of flammable anesthetic mixture with air or with oxygen or nitrous oxide other than nasal or mask type.
  - Always make sure that there are no water leaks before starting a defibrillation.
  - Always monitor the patient for shivering, temperature, signs of intolerance, and skin condition when using this product.
- 

## Starting Automatic therapy mode

In Automatic mode, the therapy cools or warms the patient to a selected patient target temperature. The product in automatic mode continually measures the patient temperature and automatically adjusts the water temperature until the selected patient target temperature is achieved. After the selected patient target temperature is achieved, the product maintains this temperature for the duration of the therapy.

To start Automatic therapy mode:

1. Prepare the thermal transfer devices for therapy.
2. See *Filling a thermal transfer device* (page 26).
3. Apply the thermal transfer device to the patient.
4. Connect the reusable adapter cable to port A or port B on the product. Make sure that the probe is fully seated.
5. Apply the sensing end of a patient probe to the patient based on your hospital protocol and secure the product to reduce the risk of accidental dislodgment.
6. Connect the patient temperature probe to the reusable adapter cable. See *Selecting and connecting a temperature probe* (page 20).
7. Tap to Confirm the current patient temperature.
8. Tap the Automatic therapy mode button.
9. Select the target patient temperature.
10. See *Setting or editing the cooling rates* (page 28) or *Setting or editing the warming rates* (page 28)

### Note

- The controller determines **Warming** or **Cooling** therapy based on the selected water target temperature and the current water temperature.
- Do not place additional heat sources between the patient and thermal transfer device.
- After the patient target temperature is achieved, patient temperature is controlled within  $\pm 0.3^{\circ}\text{C}$ .
- If the patient temperature is not within  $0.5^{\circ}\text{C}$  of the current target temperature, the yellow patient icon will flash and the patient temperature deviation alarm will sound. This occurs after the initial patient target temperature is achieved.

## Setting or editing the cooling rates

Setting cooling rates is for Automatic Mode only.

1. To set the cooling temperature, highlight your choice of cooling rates.



| Select a cooling rate                                                                 | Description                                                                                                                       |
|---------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------|
|  Max | Maximum: approaches patient target temperature as fast as possible                                                                |
|  Med | Medium: water is cooled to target, with a max of $15.0^{\circ}\text{C}$ difference between the patient and the water temperature  |
|  Min | Minimum: water is cooled to target, with a max of $10.0^{\circ}\text{C}$ difference between the patient and the water temperature |

2. Tap Confirm.
3. Tap the Edit button to make changes.

## Setting or editing the warming rates

Setting warming rates is for Automatic Mode only.

1. To set the warming temperature, highlight your choice of warming rates.



| Select a warming rate                                                                          | Description                                                                                                                                                                                       |
|------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  Max        | Maximum: approaches patient target temperature as fast as possible                                                                                                                                |
|  Med        | Medium: patient temperature increases at a rate of $4.0^{\circ}\text{C}$ in 12 hours ( $0.33^{\circ}\text{C/hour}$ ).                                                                             |
|  Min        | Minimum: patient temperature increases at a rate of $4.0^{\circ}\text{C}$ in 24 hours ( $0.17^{\circ}\text{C/hour}$ ).                                                                            |
|  Set Custom | Set Custom: patient temperature increases at a customized temperature and time period the operator selects. The temperature increases $0.05^{\circ}\text{C/hour}$ to $0.5^{\circ}\text{C/hour}$ . |

2. If you select Set Custom, tap the Increase and Decrease buttons to set the rate (Figure 19).




Figure 19 – Set custom warming rate

3. Tap Confirm.
4. Tap the Edit button to make changes.

## Starting Manual mode

In Manual mode, the therapy will cool or warm the water to a selected water target temperature. The operator must observe the patient's temperature and manually adjust the water temperature to obtain the desired patient temperature.

1. If desired, select and place the sensing end of the patient probe based on hospital protocol. Connect the reusable adapter cable to port A or port B on the product. See *Selecting and connecting a temperature probe* (page 20).
2. Prepare the thermal transfer devices to be used for therapy.
3. See *Filling a thermal transfer device* (page 26).
4. Apply the thermal transfer device to the patient.
5. Tap Manual mode. The default water target temperature is 40.0° C (104.0° F) upon initial entry.
6. Tap Confirm.
7. To select the desired water temperature, tap the Increase or Decrease buttons or hold the button to go faster.
  - a. To edit water temperature, tap the Edit button.
8. Tap Confirm.

### Note

- The controller determines **Warming** or **Cooling** therapy based on the selected water target temperature and the current water temperature.
- In Manual mode, only the water temperature is controlled.
- A temperature probe is not required when operating in Manual mode.
- After the water target temperature is achieved, water temperature is controlled within +/-0.3° C.

## Starting Monitor mode

In Monitor mode, no therapy will be delivered, only a display of the current patient temperature.

To start Monitor mode:

1. Connect the reusable adapter cable to port A or port B on the controller. Make sure that the probe is fully seated.
2. Apply the sensing end of patient probe to patient based on hospital protocol. Secure the patient probe to reduce the risk of accidental dislodgment.
3. Tap the Monitor button.
4. Connect the patient temperature probe to the end of the reusable adapter cable. See *Selecting and connecting a temperature probe* (page 20).

**Note** - If the product senses a patient probe temperature below 36.0° C (96.8° F) or above 38.0° C (104° F), the normothermia alarm displays and an audible alarm sounds.

5. Tap Confirm. The screen will display the current patient temperature.

## Switching modes

Tap Edit and select a different therapy mode.

## Pausing and resuming therapy

To pause therapy, press and hold the Pause Therapy button for two seconds.



To resume therapy, press and hold the Pause Therapy button for two seconds.


## Displaying the data storage

The system gathers data at 60 second intervals and is limited to 90 minutes of storage. The graphical display defaults to show data for all four variables in the Manual and Automatic modes.

To display the patient data graph:



1. Tap the Graph icon.



**Figure 20 – Graphical display**

- Primary patient temperature reading from attached probe (A) (Figure 20)
  - Intermediate target temperature (B)
  - Water temperature (C)
  - Power level (D)
2. To view data values, hide values, or data lines, tap an icon until the data you desire appears for the selected icon.
  3. Tap next to see the current values for each variable.
  4. To exit, tap the Graph icon or the Exit button.

### Note

- In monitor mode, only the patient temperature data is displayed (A).
- The graph icon is only available when a therapy is active.
- The patient data remains until the product sleeps or you powered off the product.
- If you lose power for more than 10 minutes, the data will not be retrievable.

## Opening and securing items in the storage compartment

The storage compartment holds a maximum of 3 lb (1.36 kg).

To open the storage compartment door, lift up on the door (Figure 21).

The storage compartment secures the following items:

- Two patient probes
- Two reusable adapter cables
- One reusable patient temperature output cable
- Product operations manual




Figure 21 – Storage compartment

### Note

- Make sure that the items are securely inside the compartment and not blocking the magnets.
- When closing the storage compartment door, do not place your fingers between the storage compartment door and the sides of the controller.

## Stopping therapy or powering off the product

To stop therapy or power off the controller:

1. Press and hold the Stand-by button for two seconds.
2. Unplug the product from the wall outlet.

**Note** - If storing the product, see *Storing the controller* (page 34).

## Draining the thermal transfer devices

Read the manufacturer's operations manual for the individual thermal transfer devices (blankets and wraps) for warnings, cautions, and safe operating instructions before use. Make sure that you drain the hoses before you put them in storage.

1. Unplug the power cord from the wall outlet.
2. Remove the thermal transfer device from the patient.
3. Open the clamps on the hoses and thermal transfer devices, if applicable. See Figure 15.
4. Raise the thermal transfer devices attached to the hose above the ports on the controller. Gravity helps to drain the water into the controller.
5. Allow most of the water to drain back into the controller. (Approximately 10 minutes).
6. See *Connecting and disconnecting thermal transfer devices* (page 22).
7. See *Disconnecting the insulated hoses* (page 22).
8. See *Storing the power cord and hoses* (page 34).

9. Discard the disposable thermal transfer devices based on your local waste management protocol.

## Draining water from the reservoir

To drain the water from the reservoir:


1. See *Removing and replacing the reservoir* (page 24).
2. Dispose of the water per hospital protocol.
3. Replace the reservoir.

**Note** - Make sure that the reservoir is dry before you store the product.

## Draining water from the controller and hoses


Make sure that the controller and all components are dry before you store the product. Make sure to drain the hoses before you store them.

1. Place the controller over a floor drain or drain pan.
2. Remove the reservoir and pull up on the controller drain plug (A) to open the drain (Figure 22).



**Figure 22 – Drain plug**

3. Connect a hose to each port.
  - a. If you have Colder style connector hoses, attach the service tool adapter hose (8001-999-017).
  - b. If you have Klik-Tite hoses, make sure that the connectors and clamps are open (Figure 23).



**Figure 23 – Klik-Tite open**

4. Raise all the hoses completely above the connection ports on the controller.
5. Allow the product to drain for a minimum of two minutes.

6. Push down on the drain plug to close the drain.
7. Empty the water from the drain pan, if applicable.
8. Replace the reservoir.

## Emptying water with a drain pan

To empty water with a drain pan:

1. Slide an empty drain pan under the front of the controller.




Figure 24 – Slide under the controller




Figure 25 – Drain pan in place

2. Allow the controller to drain for two minutes.




Figure 26 – Drain for two minutes




Figure 27 – Remove to drain and store

3. Remove the drain pan slowly.

4. Empty the water from the drain pan into the reservoir for disposal.
5. Dry the drain pan before you store for future use.

## Storing the power cord and hoses

After you complete therapy or when you transport the product, store the power cord and hoses.

---

### CAUTION

- Do not hang items on the controller handle to avoid the risk of tipping the product.
  - Always store the power cord, cables, and hoses before you transport the product to reduce the risk of a trip hazard.
- 

To store the power cord and hoses:

1. Connect the ends of the connector hoses together, if applicable.
2. Coil and fasten the hose with the management straps (Figure 28).
3. Unplug the power cord from the wall outlet.
4. Store the power cord with the management straps (Figure 28).



Figure 28 – Management straps

## Storing the controller

Storage is equal to or greater than 7 days without use.

---

### CAUTION

- Do not store the product with water in the device.
  - Always store the product within the specified environmental condition values.
- 

To store the controller:

1. See *Disinfecting the internal water circuit and hoses every 14 days* (page 43).
2. See *Draining the thermal transfer devices* (page 31).
3. See *Cleaning the external surfaces* (page 40).
4. See *Disinfecting external surfaces* (page 41).

### Note

- Always bring the controller to room temperature after high or low temperature storage.
- Always store the controller with the reservoir in place.

## Transporting the product


Make sure to follow these procedures for transporting the product to avoid the risk of possible injury or equipment damage.

---

### CAUTION

- Always use extra care when you transport the product long distances and on inclines greater than five degrees. Ask for help, if necessary, to avoid the risk of tipping.
  - Always use the handle to move the product. Do not attempt to move the product by pulling on cables, hoses, or by any other means.
  - Avoid ramps that are steeper than ten degrees to avoid tipping the product.
  - Do not hang items on the controller handle to avoid the risk of tipping the product.
  - Always store the power cord, cables, and hoses before you transport the product to reduce the risk of a trip hazard.
- 

1. Make sure that the pathway is clear.
2. Unplug the power cord from the wall outlet. See *Storing the power cord and hoses* (page 34).
3. Make sure to place the product with the ports facing the front (Figure 29).



**Figure 29 – Transport position**

4. Use the handle to push the product.
5. Limit movement to a slow careful walk.
6. Use two or more operators to move the system on inclines or long distances.

### Note

- Wheel chair ramps are usually less than five degrees.
- The system weighs 150 lb (68.0 kg) when dry. Weight also depends on other items in the storage compartment.

## Accessories and parts

These accessories may be available for use with your product. Confirm availability for your configuration or region. Call 3T Medical Customer Service: (734) 525-9580. For more information, see the thermal transfer device instructions for use.

### Thermal transfer devices

| Thermal transfer device                                                        | Connector type | Part number  | Size                 |                  |
|--------------------------------------------------------------------------------|----------------|--------------|----------------------|------------------|
| Rapr-Round, small/medium vest chest                                            | Clik-Tite      | 8001-061-530 | 32 in. to 46 in.     | 81 cm to 117 cm  |
| Rapr-Round, large vest chest                                                   | Clik-Tite      | 8001-061-535 | 46 in. to 54 in.     | 117 cm to 137 cm |
| Rapr-Round, leg wrap, one size for the left or right leg - thigh circumference | Clik-Tite      | 8001-061-540 | 20.5 in. to 28.5 in. | 52 cm to 72 cm   |
| Mul-T-Blanket, adult                                                           | Colder         | 8001-061-610 | 25 in. x 64 in.      | 64 cm x 163 cm   |
| Mul-T-Blanket, pediatric                                                       | Colder         | 8001-061-612 | 22 in. x 33 in.      | 56 cm x 84 cm    |
| Rapr-Round, small/medium vest chest                                            | Colder         | 8001-061-630 | 32 in. to 46 in.     | 81 cm to 117 cm  |
| Rapr-Round, large vest chest                                                   | Colder         | 8001-061-635 | 46 in. to 54 in.     | 117 cm to 137 cm |
| Rapr-Round, leg wrap, one size for the left or right leg - thigh circumference | Colder         | 8001-061-640 | 20.5 in. to 28.5 in. | 52 cm to 72 cm   |
| Mul-T-Blanket, adult                                                           | Clik-Tite      | 8001-061-810 | 25 in. x 64 in.      | 64 cm x 163 cm   |
| Mul-T-Blanket, pediatric                                                       | Clik-Tite      | 8001-061-812 | 22 in. x 33 in.      | 56 cm x 84 cm    |
| Altrix Temperature Management Wraps                                            |                |              |                      |                  |
| Small torso                                                                    | Colder         | 8003-002-001 | 26 in. - 47 in.      | 66 cm - 119 cm   |
| Large torso                                                                    | Colder         | 8003-002-002 | 35 in. - 60 in.      | 89 cm - 152 cm   |
| Small left thigh                                                               | Colder         | 8003-002-003 | 16 in. - 28 in.      | 41 cm - 71 cm    |
| Small right thigh                                                              | Colder         | 8003-002-004 | 16 in. - 28 in.      | 41 cm - 71 cm    |
| Large left thigh                                                               | Colder         | 8003-002-005 | 19.5 in - 32 in.     | 50 cm - 81 cm    |
| Large right thigh                                                              | Colder         | 8003-002-006 | 19.5 in. - 32 in.    | 50 cm - 81 cm    |

### Thermal transfer device kits

| Kit part number | Contents     | Quantity | Type of Connector |
|-----------------|--------------|----------|-------------------|
| 8001-061-550    | 8001-061-530 | 1        | Clik-Tite         |
|                 | 8001-061-540 | 2        |                   |
| 8001-061-560    | 8001-061-535 | 1        |                   |
|                 | 8001-061-540 | 2        |                   |

| Kit part number | Contents     | Quantity | Type of Connector |
|-----------------|--------------|----------|-------------------|
| 8001-061-650    | 8001-061-630 | 1        | Colder            |
|                 | 8001-061-640 | 2        |                   |
| 8001-061-660    | 8001-061-635 | 1        |                   |
|                 | 8001-061-640 | 2        |                   |
| 8003-003-001    | 8003-002-003 | 1        | Colder            |
|                 | 8003-002-004 | 1        |                   |
|                 | 8003-002-001 | 1        |                   |
| 8003-003-002    | 8003-002-005 | 1        | Colder            |
|                 | 8003-002-006 | 1        |                   |
|                 | 8003-002-002 | 1        |                   |
| 8003-003-003    | 8003-002-003 | 1        | Colder            |
|                 | 8003-002-004 | 1        |                   |
| 8003-003-004    | 8003-002-005 | 1        | Colder            |
|                 | 8003-002-006 | 1        |                   |

## Patient temperature probes

| Patient temperature probes                     | Part Number  |
|------------------------------------------------|--------------|
| Adhesive skin temperature sensing probe        | 8001-063-401 |
| 9FR General purpose temperature sensing probe  | 8001-063-409 |
| 12FR General purpose temperature sensing probe | 8001-063-412 |
| 14FR Foley catheter temperature sensing probe  | 8001-063-414 |
| 16FR Foley catheter temperature sensing probe  | 8001-063-416 |

## Cables

| Description                               | Part Number  |
|-------------------------------------------|--------------|
| Reusable adapter cable                    | 8001-064-110 |
| Reusable patient temperature output cable | 8001-064-120 |

## Hoses

| Description                        | Part Number  |
|------------------------------------|--------------|
| Insulated Klik-Tite Hose           | 8001-064-035 |
| Insulated Colder Connector Hose    | 8001-064-135 |
| Altrix Temperature Management Hose | 8003-002-007 |

## Cleaning tools

| Description               | Part Number  |
|---------------------------|--------------|
| Klorkleen Medical tablets | 8001-999-228 |
| Service tool adapter hose | 8001-999-017 |
| Drain pan                 | 800109990234 |

# Preventive maintenance

Remove the product from service before you perform the preventive maintenance inspection. Check all items listed during annual preventive maintenance for all 3T Medical products. You may need to perform preventive maintenance checks more often based on your level of product usage. Service only by qualified service technician.

**Note** - Complete functional checks as indicated in the maintenance manual.

Inspect all of the following items:

- \_\_\_\_\_ All fasteners secure
- \_\_\_\_\_ Power cord and plug for fraying
- \_\_\_\_\_ Condition of covers and push handle for damage
- \_\_\_\_\_ Hose ports are operational
- \_\_\_\_\_ Ground chain attached
- \_\_\_\_\_ LCD is not cracked
- \_\_\_\_\_ LCD operational
- \_\_\_\_\_ Touch screen operational
- \_\_\_\_\_ Wheels for smooth operation
- \_\_\_\_\_ Rear caster wheels for free swivel action
- \_\_\_\_\_ Both rear wheels lock secure when applied
- \_\_\_\_\_ Front and rear wheels are not loose or wobbly
- \_\_\_\_\_ Battery backup, functional check
- \_\_\_\_\_ Alarm system - visual and audible, functional check
- \_\_\_\_\_ High temp cutout checks
- \_\_\_\_\_ Water temperature and flow verification, functional check
- \_\_\_\_\_ Probe resistance, functional check
- \_\_\_\_\_ Clear RFU codes
- \_\_\_\_\_ Ground impedance not more than 100mΩ (milliohms)
- \_\_\_\_\_ Current leakage not more than 300 (microamps)
- \_\_\_\_\_ Integrity of all clamps and clamped joints located in the air elimination circuit

Replace the following on an annual basis:

- \_\_\_\_\_ 9 V battery
- \_\_\_\_\_ Condenser inlet filter
- \_\_\_\_\_ Air eliminator hose
- \_\_\_\_\_ Collars
- \_\_\_\_\_ Inlet filter

|                        |
|------------------------|
| Product Serial Number: |
| Completed by:          |
| Date:                  |

# Cleaning

## Cleaning the external surfaces

Clean the external surfaces of the controller and system components before each use. System components may be subject to contamination during use from contact with soiled hands of the user, airborne pathogens, and unexpected or accidental events. Make sure you remove all visible soils.

---

**CAUTION** - Do not power wash this product.

---

### Tools Required:

- Mild soap
- Soft, lint free cloth (2 or more)

Validated mild soap:

- Enzol® Enzymatic Instrument Cleaner by Johnson & Johnson

To clean the external surfaces of the controller and system components:

See *Product illustration* (page 9) for clarification of product component names and locations.

1. Unplug the power cord from the wall outlet.
2. Apply wheel locks.
3. Undo power cord and hose straps.
4. Unravel and lay out hoses, cables, and power cord.
5. Disconnect the hoses. Push back on the retaining collar of the port on the controller. Pull the hose to disconnect.
6. Disconnect the patient temperature probe cable from the port.
7. Remove the reservoir. Pull forward at an angle, and lift out the reservoir.
8. If necessary, empty the water from the reservoir. Dispose of the water per your hospital protocol.
9. Prepare a mild soap and water solution as described by the manufacture.
10. Wipe the inside and outside of the reservoir and reservoir lid, with a soft, lint free cloth moistened with soap and water solution.
11. Wipe the hoses and patient temperature probe cables, with a soft, lint free cloth moistened with soap and water solution.
12. Wipe the controller surfaces while the reservoir is removed with a soft, lint free cloth moistened with soap and water solution. Also wipe the following system components:
  - Hose connectors
  - Power cord
  - Hose and power cord management straps
  - Storage compartment door
  - Inside storage compartment
  - Graphical user interface display
  - Handle
13. Wipe the controller, reservoir and reservoir lid surfaces, and system components with a clean, dry cloth to remove any excess liquid.
14. Replace the reservoir.
15. Allow the external surfaces of the controller and components to dry thoroughly.

# Disinfecting

## Disinfecting external surfaces

Disinfect the external surfaces of the controller and system components before each use. System components may be subject to contamination during use from contact with soiled hands of the user, airborne pathogens, and unexpected or accidental events. Follow your hospital protocols for disinfecting the product. Make sure to follow the manufacturer's instructions for the disinfectants.

---

**CAUTION** - Do not use quaternaries that contain glycol ethers as they may damage the reusable accessories.

---

**Note** - If the product is visibly soiled, clean the surface before disinfecting.

### Tools Required:

- Personal protection equipment (PPE) as recommended by the disinfectant manufacturer's instructions
- Soft, lint free cloth (2 or more)
- Disinfectant
- 2 gallons (7.6 L) of sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns

Recommended disinfectants for the external surface of the controller and system components:

- Quaternary cleaners (active ingredient - ammonium chloride)
- Phenolic cleaners (active ingredient - o-phenylphenol)
- Chlorinated bleach solution ((1 part bleach solution (5.25% sodium hypochlorite) to 100 parts of water which equals 520 ppm available chlorine (40 ml of a 5.25% bleach solution per 4000 ml water))

Validated disinfectants for the external surface of the controller and system components:

- Sodium hypochlorite based - **Clorox®** Healthcare Bleach Germicidal Cleaner (EPA registration number 56392-7)

To disinfect the external surfaces of the controller and system components:

See *Product illustration* (page 9) for clarification of product component names and locations.

1. Use PPE as recommended by the disinfectant manufacturer's instructions.
2. Unplug the power cord from the wall outlet.
3. Apply the wheel locks.
4. Unfasten the power cord and hose straps.
5. Unravel and lay out hoses, cables, and power cord.
6. Disconnect the hoses. Push back on the retaining collar of the port on the controller. Pull the hose to disconnect.
7. Disconnect the patient temperature probe cable.
8. Remove the reservoir. Pull forward at an angle, and lift out the reservoir.
9. If necessary, empty the water from the reservoir. Dispose of the water per your hospital protocol.
10. Prepare disinfectant solution as described by the manufacture.
11. Apply disinfectant solution to the inside and outside of the reservoir and reservoir lid, with a soft, lint free cloth moistened with disinfectant. Reapply disinfectant to cloth as needed.
12. Apply disinfectant solution to the hoses and patient temperature probe cables, with a soft, lint free cloth moistened with disinfectant. Reapply disinfectant to cloth as needed.
13. Apply disinfectant solution to the controller surfaces while the reservoir is removed with a soft, lint free cloth moistened with disinfectant. Reapply disinfectant to the cloth as needed. Also wipe the following system components:
  - Hose connectors

- Power cord
- Hose and power cord management straps
- Storage compartment door
- Inside storage compartment
- Graphical user interface display
- Handle

14. Follow specified contact time in accordance with the disinfectant manufacturer's instructions for use.
15. To rinse, wipe the hoses and patient temperature probe cables, with a soft, lint free cloth moistened with sterile distilled water.
16. To rinse, wipe the controller, reservoir and reservoir lid surfaces, and system components with a lint free cloth moistened with sterile distilled water.
17. To dry, wipe the controller, reservoir, reservoir lid surfaces, and system components with a clean, dry cloth to remove any excess liquid.
18. Replace the reservoir.
19. Allow the external surfaces of the controller and components to dry thoroughly.
20. Store the power cord, cables, and hoses.

## Disinfecting the internal water circuit and hoses every 14 days

Use the **Klorkleen Medical** disinfectant tablets manufactured by **Medentech** (EPA registration number 71847-2-106) before first use, at least every 14 days, and before storage. **Klorkleen Medical** disinfectant tablets have been validated for internal water circuit disinfection. Make sure that you follow the disinfectant manufacturer's guidelines to avoid the risk of injury. Failure to follow the disinfectant's instructions may void your warranty.

---

### CAUTION

- Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.
  - Do not disinfect the internal water system with a thermal transfer device attached as this may cause a leak.
  - Do not use bleach or any other cleaning or disinfectant agents for internal circuits. This could result in damage to the product. Only use approved disinfectant tablets.
  - Always drain the product before disinfecting the internal water circuit. Failure to drain the product may reduce the effectiveness of the disinfection process.
- 

**Note** - Disinfection of the Altrix internal water system was validated using *M. mucogenicum*.

### Tools Required:

- 2 gallons (7.6 L) of sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns
- Personal protection equipment (PPE) as recommended by the disinfectant manufacturer's instructions
- Soft, lint free cloth (2 or more)
- 2 **Klorkleen Medical** disinfectant tablets (active ingredient is sodium dichloroisocyanurate (NaDCC); 1874 ppm disinfectant solution when dissolved in one gallon of water)
- Service tool adapter hose (8001-999-017) for **Colder** style connector hoses
- Floor drain or drain pan (800109990234)


See *Product illustration* (page 9) for clarification of product component names and locations.

## Draining the internal water circuit and hoses for disinfection

1. Unplug the power cord from the wall outlet.
2. Place the controller over a floor drain or drain pan.

**Note** - For best results, the floor drain should be within reach of a wall outlet to power on the controller.

3. To drain the controller, pull up on the controller drain plug (A) to open the drain (Figure 30). Leave the drain open.



**Figure 30 – Drain plug**

4. Connect a hose to each port (Figure 31).




**Figure 31 – Hoses connected**

5. Close the connector ends of all three hoses:
  - a. If you have **Colder** style connector hoses, attach the service tool adapter hose (8001-999-017) (Figure 32). Complete this for all three hoses.



**Figure 32 – Colder style connector hose connected to a tool adapter hose**

- b. If you have **Clik-Tite** hoses, make sure that the connector ends are connected and closed (A), and clamps are open (B). Complete this for all three hoses. Figure 33



**Figure 33 – Clik-Tite hose ends are closed and clamps are open**

6. To fully drain the hoses, raise all the hoses (Figure 34) above the connection ports on the controller.

**Note** - For best performance, hang the hoses to keep them raised. Do not lower the hoses until you have completed the disinfection and rinsing process.



**Figure 34 – Raise the hoses**

7. Allow the controller and hoses to drain for a minimum of two minutes.
8. Push down on the drain plug to close the drain.
9. Empty the water from the drain pain, if applicable.

## Disinfecting the internal water circuit and hoses

**CAUTION** - Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.

1. Use personal protection equipment as recommended by the **Klorleem Medical** disinfectant manufacturer's instructions for use.
2. Put two **Klorleem Medical** disinfectant tablets into the reservoir.
3. Using appropriate measuring equipment, fill the empty reservoir with one gallon (3.8 L) of sterile distilled water.

**Note** - Always allow the disinfectant tablets to completely dissolve before starting the 20 minute disinfection cycle.

4. Place the reservoir into the controller.
5. Disconnect the bottom hose from the bottom right port (Figure 35).



**Figure 35 – Disconnected hose**

6. Connect the bottom hose end to the hydraulic connector in the lid of the reservoir (Figure 36).



**Figure 36 – Bottom hose end in the lid of the reservoir**

7. Plug the power cord into a wall outlet.

8. Press and hold the  Stand-by button.

9. Tap the Manual mode icon. 

10. Tap Confirm. 

11. Set the water target temperature to 25.0° C (77.0° F).

12. Tap Confirm. 

13. Allow the controller to run for 2 minutes to purge air from the system.

14. Run the controller for 20 minutes.  00:22:00

15. After 22 minutes, turn the controller off by pressing and holding the Stand-by button for two seconds. 

16. Unplug the power cord from the wall outlet.

17. Place the controller over a floor drain.


18. Remove the reservoir. Pull forward at an angle, and lift out the reservoir.

19. Remove the bottom hose end from the hydraulic connector adapter in the reservoir lid by pushing down on the collar.

20. Empty water from the reservoir, dispose of the water per hospital protocol.

**Note** - Do not rinse the reservoir.

21. Pull up on the controller drain plug (Figure 37) to open the drain.



**Figure 37 – Drain plug**

22. Make sure that all three hoses remain raised above the connection ports for draining.

23. Allow the controller and hoses to drain for a minimum of two minutes.

24. Push down on the controller drain plug to close the drain.

25. When the controller and hoses are drained, continue to *Rinsing the internal water circuit and hoses* (page 47).

## Rinsing the internal water circuit and hoses

---

**CAUTION** - Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.

---

1. Using appropriate measuring equipment, fill the empty reservoir with one gallon (3.8 L) of sterile distilled water.
2. Place the reservoir into the controller.
3. Connect the bottom hose end to the hydraulic connector in the lid of the reservoir.



**Figure 38 – Bottom hose end in the lid of the reservoir**

4. Plug the power cord into a wall outlet.

5. Press and hold the  Stand-by button.

6. Tap the Manual mode icon. 

7. Tap Confirm. 

8. Select the water target temperature of 25.0° C (77.0° F).

9. Tap Confirm. 

10. Allow the controller to run for 2 minutes to purge air from the system.

11. Allow the controller to run for 5 minutes.  00:07:00

**Note** - The timer will run on the main display, follow the current therapy duration timer.

12. After seven minutes, turn the controller off by pressing and holding the Stand-by button for two seconds. 
13. Unplug the power cord from the wall outlet.
14. Place the controller over a floor drain or drain pan.
15. Remove the reservoir. Pull forward at an angle, and lift out the reservoir.
16. Remove the bottom hose end from the hydraulic connector adapter in the reservoir lid by pushing down on the collar.
17. Empty water from the reservoir, dispose of the water per hospital protocol.
18. Pull up on the controller drain plug to open the drain.
19. Make sure that all three hoses remain raised above the connection ports for draining.
20. Allow the controller and hoses to drain for a minimum of two minutes.
21. Push down on the controller drain plug to close the drain.
22. Empty the water from the drain pan, if applicable.
23. Wipe the inside and outside of the reservoir and reservoir lid, with a dry, soft, lint free cloth.
24. Place the reservoir into the controller.
25. Disconnect and store the service tool adapter hoses from all three of the hoses. (If applicable, when used with **Colder** style hoses.)
26. Store the power cord, cables, and hoses.

# Troubleshooting

**CAUTION** - Always use sterile distilled water or distilled water that has been passed through a filter less than or equal to 0.22 microns with this product.

| Problem                                                                    | Possible cause                                                                | Action                                                                                                                                                                                                                                    | Remove from use (RFU)                                           |
|----------------------------------------------------------------------------|-------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------|
| Controller will not turn on                                                | Power cord is not plugged into a properly grounded hospital grade wall outlet | Insert the plug fully into the properly grounded hospital grade wall outlet                                                                                                                                                               | If the product does not turn on after trying a different outlet |
|                                                                            | Damaged power cord or plug                                                    | Visually make sure that the power cord is not damaged                                                                                                                                                                                     | If damaged, RFU                                                 |
| Controller user interface blackout                                         | Power outage                                                                  | If the Stand-by button is solid green, visually inspect the LCD for damage                                                                                                                                                                | If damaged, RFU                                                 |
| Product alarming, user interface blackout                                  | Power outage                                                                  | If the Stand-by button is yellow and flashing, visually inspect the LCD for damage                                                                                                                                                        | If damaged, RFU                                                 |
| Temperature probe                                                          | Not responding, does not connect, temperature outside of range                | Replace temperature probe<br>Check connections                                                                                                                                                                                            | If damaged, RFU                                                 |
| Controller will not heat                                                   | Reservoir is empty                                                            | Fill reservoir, tap Confirm that sterile distilled water has been added, restart                                                                                                                                                          | If filling the reservoir does not resolve                       |
| Controller will not cool                                                   | Reservoir is empty                                                            | Fill reservoir, tap Confirm that sterile distilled water has been added, restart                                                                                                                                                          | If filling the reservoir does not resolve                       |
| Thermal transfer device not filling with water or ports not detecting flow | Locking ring on Clik-Tite connector is not snapped into place                 | Check the <b>Clik-Tite</b> connection<br><br>Replace the cable or thermal transfer device<br>Thermal transfer device may be too high, lower the bed level<br>Thermal transfer device may be folded, lay flat to make sure the water flows | Not applicable                                                  |
|                                                                            | Quick disconnect not seated properly                                          | Secure the thermal transfer device connection to the controller<br><br>Replace the cable or thermal transfer device                                                                                                                       | Not applicable                                                  |
| Water level alarm                                                          | Water level too low                                                           | Fill reservoir                                                                                                                                                                                                                            | Not applicable                                                  |

| Problem                          | Possible cause                                                                                                                                                                                                                         | Action                                                                                                                                                                                                                                                                                                                                                                                                                      | Remove from use (RFU) |
|----------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|
| Patient temperature              | Out of range                                                                                                                                                                                                                           | Check the placement of the probe                                                                                                                                                                                                                                                                                                                                                                                            | Not applicable        |
| Patient temperature output (PTO) | External device output displays high value >45° C when input is out of range (As a result of one of the following: patient probe disconnected, controller in standby/sleep mode, patient temperature outside the range of 25° - 45°C). | To resume calibration: <ul style="list-style-type: none"> <li>• Disconnect the external device from the reusable adapter cable</li> <li>• Tap the Help button to display the alarm screens</li> <li>• Find the "Temperature Output Alarm" screen</li> <li>• Tap Confirm to restart calibration</li> <li>• Wait until the "Monitor" is solid</li> <li>• Connect the external device to the reusable adapter cable</li> </ul> | Not applicable        |

# Alarm conditions

The alarm rank establishes the order of presentation of the alarm on-screen message. The D in the table indicates the alarm is deactivated during that mode. Maintenance and RFU modes are always in the deactivated mode and are not listed in the table.

This product maintains the individual alarm status for all alarms as defined below.

- Alarm condition present
- Visual indicator state
- On screen message
- Audible indicator state
- Current timer for audio pause activation
- Alarm rank per therapy mode

| Alarm (med)                     | Stand-by | Auto | Auto Paused | Manual | Manual Paused | Monitor |
|---------------------------------|----------|------|-------------|--------|---------------|---------|
| Remove from use                 | 0        | 0    | 0           | 0      | 0             | 0       |
| Power loss                      | D        | 1    | 1           | 1      | 1             | 1       |
| Check patient probe             | D        | 7    | 7           | 11     | 11            | 4       |
| Patient probe malfunction       | D        | 6    | 6           | 10     | 10            | 3       |
| Probe disconnected              | D        | 5    | 5           | 9      | 9             | 2       |
| Patient temperature deviation   | D        | 9    | 9           | D      | D             | D       |
| Water temperature deviation     | D        | D    | D           | 7      | D             | D       |
| Check water flow (all ports)    | D        | D    | D           | 6      | D             | D       |
| No water flow alarm (all ports) | D        | 4    | D           | 4      | D             | D       |
| No water                        | D        | 2    | D           | 2      | D             | D       |
| Therapy paused timed out        | D        | D    | 3           | D      | 3             | D       |
| Alarm (low)                     | Stand-by | Auto | Auto Paused | Manual | Manual Paused | Monitor |
| Normothermia deviation          | D        | D    | D           | D      | D             | 5       |
| Patient temperature deviation   | D        | 12   | 12          | D      | D             | D       |
| Power backup level              | 1        | 19   | 19          | 14     | 14            | 6       |
| Patient output deviation        | D        | 22   | 22          | 17     | 17            | 7       |

## Note

- If more than one alarm is active at the same time, the product maintains the active state for the individual alarm including the audio pause timer. The screen alarms are displayed with the highest level alarm first with a page toggle to allow the operator to scroll to the subsequent alarms.
- The Paused in Auto Paused and Manual Paused refer to the Therapy Pause state.

## Check patient probe alarm

This alarm notifies the operator that data provided by the probe is not normal or appears removed.

The product only activates the Check Patient Probe Alarm when met during an Active Therapy. Otherwise, the alarm is disabled.

Alarm generation:

Primary patient temperature changes by more than 1.0° C within two minutes.

### Note

- The product will deactivate the heat exchange and keep the pump activated as requested by the Active Therapy.
- The audio paused timeout is 5 minutes for the check patient probe alarm.

## Patient probe malfunction alarm

This alarm notifies the operator that the probe is not providing information to the product during an active therapy.

Alarm generation:

When the primary patient probe is in a shorted, opened condition, or out of range for more than 30 seconds, the product will display the Patient Probe Malfunction Alarm.

### Note

- The product will deactivate the heat exchange and keep the pump activated as requested by the Active Therapy.
- The audio paused timeout is 5 minutes for the patient probe malfunction alarm.

## Patient probe disconnect alarm

This alarm notifies the operator that the probe is not providing information to the product during an active therapy.

Alarm generation:

When the adapter cable for the primary probe is removed and the reading of the Primary Patient probe is out of range for more than 30 seconds, the product displays the Patient Probe Disconnected alarm.

**Note** - The audio paused timeout is 5 minutes for the patient probe disconnected alarm.

## Patient temperature deviation medium alarm

This alarm notifies the operator that the patient is not responding as expected in the active therapy.

Alarm generation:

The product will display the Patient Temperature Deviation Medium Alarm, if after the initial achievement of the current patient temperature target, the actual primary patient temperature becomes more than 1.0° C above or below the Current Target Temperature.

**Note** - The audio paused timeout is 30 minutes for the patient temperature deviation medium alarm.

## Patient temperature deviation low alarm

This alarm notifies the operator that the patient is not responding as expected in the active therapy.

Alarm generation:

The product will display the Patient Temperature Deviation Low Alarm, if after the initial achievement of the current patient temperature target, the actual primary patient temperature becomes more than 0.5° C above or below the Current Target Temperature.

**Note** - The low priority audible alarm is repeated every 30 minutes.

## Patient temperature output deviation alarm

This alarm notifies the operator that the patient temperature output is out of range or there is a calibration error.

Alarm generation:

The product will display the Patient Temperature Output Deviation Alarm, if the calibration has failed, or the patient temperature output is out of range.

**Note** - The audio paused timeout is 5 minutes for the no water flow alarm.

## Normothermia deviation alarm

This alarm notifies the operator of that the Primary Patient Temperature is out of range.

Alarm generation:

If the actual Primary Patient Temperature is lower or equal to 35.9° C or higher or equal to 38.1° C, the controller will display the Patient Normothermia Deviation alarm.

**Note** - The low priority audible alarm is repeated every 30 minutes.

## Water temperature deviation alarm group

This alarm notifies the operator that the water is not responding as expected to the therapy. The product is at full power, with the current mode and the temperature selection. The water temperature cannot remain within a range of  $\pm 0.8^{\circ}\text{C}$  of the selected Water Target Temperature.

Alarm generation:

1. If the actual water temperature is more than 0.8° C above or below the Final Target Temperature for 5 minutes, the product will display the Water Temperature Deviation Alarm.
2. When the product is entering the Manual mode or you change Target Temperature, the product will pause the audible component of the Water Temperature Deviation Alarm for 4 hours. The 4 hour pause automatically cancels after the Water Temperature becomes equal to the Final Target Temperature.

**Note** - The audio paused timeout is 10 minutes for the water temperature deviation alarm.

## Check water flow alarm

This alarm notifies the operator of the quality of the flow in each individual water circuit.

Alarm generation:

- When in Manual Mode and several ports are in use for the therapy.

- You have selected an outlet port and the flow is lower than 0.8l /min for a period of 10 minutes or more. The product will display a Check Water Flow Alarm for the given port.
- Good flow is reduced to Low flow for > 10 minutes or reduced to No flow for > 15 seconds.

#### **Note**

- The alarm displays if the flow is not at an optimal level on each port. This alarm will ask the operator to confirm which ports are currently in use.
- The addition of a port does not need the operator to confirm the addition.
- The removal of a port requires the operator to confirm.
- The Check Water Flow Alarm for the given outlet port stops, if the operator confirms removal.
- The audio paused timeout is 10 minutes for the check water flow alarm.

## **No water flow alarm**

This alarm notifies the operator of water flow that is not optimal.

Alarm generation:

- When all 3 ports have a flow lower than 0.6l/min for a minimum of 20 seconds, and the pump current is greater than or equal to 1.0A and the pump is at full speed.

#### **Note**

- The pump continues to operate when the No Flow Alarm is active.
- The audio paused timeout is 5 minutes for the no water flow alarm.

## **No water alarm**

This alarm notifies the operator if there is no water in the system and the pump has shut down.

Alarm generation:

- Pump is active in the Manual or Automatic mode. The pump current is lower than 1.0 A and 1 of the ports has a low flow for more than 60 seconds.
- Pump is active in the Manual or Automatic mode. The pump current is lower than 1.0 A and all 3 port flow states are no flow for more than 60 seconds.
- The pump current is less than 0.5 A for 20 seconds.

**Note** - The audio paused timeout is 5 minutes for the no water alarm.

## **Power backup level alarm**

This alarm notifies to the operator an indicator of the status of the Backup Power Level.

#### **Note**

- The indicator will remain active until a qualified technician replaces the battery.
- There is no reduction in the usability of the product. The product remains functional and a visual alarm displays.
- The product will disable the Backup Power Level Alarm in Sleep Mode. Otherwise, the alarm is enabled.

Alarm generation:

The product will display Backup Power Level Alarm when the battery level backup power is less than 100 minutes of alarms. Once activated, the Backup Power Level Alarm will remain active until you power off the product.

**Note** - The audio paused timeout is 10 minutes for the power backup level alarm.

## Therapy paused time out alarm

This alarm converts a therapy pause into an alarm if the duration of the pause is too long.

Alarm generation:

When paused for five minutes, the product will display the Therapy Paused Time Out Alarm. After you resume the current therapy, the Therapy Paused Time Out Alarm is deactivated.

**Note** - The audio paused timeout is 10 minutes for the therapy paused time out alarm.

## Remove from use mode

The Remove From Use (RFU) is a safety mode to limit operations. A fault condition prevents the product from normal functions and requires service. The controller will stop the Active Therapy and communicate to the operator that the controller is going into RFU mode.

---

**CAUTION** - Always remove the product from use before servicing any components. Contact qualified service personnel for service.

---

Depending on the remove from use (RFU) condition, text may or may not be displayed. For example, if there is a power loss.

- Water temperature probes are out of the allowed range
- Program and data checksum error
- High thermal cutout test failed
- Backup power product replacement required
- Low or over safety temperature
- Pump over current
- Compressor power fault
- Heater power fault
- Refrigerant control valve fault
- Main DC power lost
- CAN heart beat lost
- Dual safety temperature sensors do not match readings
- Dual safety temperature sensors are out of the allowed range
- Hardware watchdog heartbeat failure

**Note** - The audio paused timeout is 10 minutes for the remove from use alarm.

# EMC Information

**CAUTION** - This equipment is not intended for use in residential environments and may not provide adequate protection to radio reception in such environments.

| Guidance and manufacturer's declaration - electromagnetic emissions                                                                                                                |                                                                                         |                                                                                                                                                                                                                 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The Altrix system is intended for use in the electromagnetic environment specified below. The customer or the user of Altrix should assure that it is used in such an environment. |                                                                                         |                                                                                                                                                                                                                 |
| Emissions test                                                                                                                                                                     | Compliance                                                                              | Electromagnetic environment                                                                                                                                                                                     |
| RF Emissions<br>CISPR 11                                                                                                                                                           | Group 1                                                                                 | The Altrix system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                          |
| RF Emissions<br>CISPR 11                                                                                                                                                           | Class A                                                                                 | The Altrix system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic Emissions<br>IEC 61000-3-2                                                                                                                                                | Class A<br>220-240V/50Hz<br>220V/60Hz<br>Does not apply to 100V<br>50/60Hz or 120V/60Hz |                                                                                                                                                                                                                 |
| Voltage Fluctuations<br>Flicker Emissions<br>IEC 61000-3-3                                                                                                                         | Complies<br>220-240V/50Hz only                                                          |                                                                                                                                                                                                                 |

**Note** - The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment, for which CISPR 11 class B is normally required, this equipment might not offer adequate protection to radio frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.

| Recommended separations distances between portable and mobile RF communications equipment and the Altrix system                                                                                                                                                                                                                                                                                                                         |                                                                |                                            |                                             |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------|--------------------------------------------|---------------------------------------------|
| The Altrix system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of Altrix can help prevent electromagnetic interferences by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Altrix system as recommended below, according to the maximum output power of the communications equipment. |                                                                |                                            |                                             |
| Rated maximum output power of transmitter<br>W                                                                                                                                                                                                                                                                                                                                                                                          | Separation distance according to frequency of transmitter<br>m |                                            |                                             |
|                                                                                                                                                                                                                                                                                                                                                                                                                                         | 150 kHz to 80 MHz<br>$D=(1.2) (\sqrt{P})$                      | 80 MHz to 800 MHz<br>$D=(0.35) (\sqrt{P})$ | 800 MHz to 2.7 GHz<br>$D=(0.70) (\sqrt{P})$ |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                    | 0.12                                                           | 0.035                                      | 0.07                                        |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                     | 0.38                                                           | 0.11                                       | 0.22                                        |

| Recommended separations distances between portable and mobile RF communications equipment and the Altrix system |     |      |     |
|-----------------------------------------------------------------------------------------------------------------|-----|------|-----|
| 1                                                                                                               | 1.2 | 0.35 | 0.7 |
| 10                                                                                                              | 3.8 | 1.1  | 2.2 |
| 100                                                                                                             | 12  | 3.5  | 7   |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

| Guidance and manufacturer's declaration - electromagnetic immunity                                                                                                                 |                                                                                                                                                                            |                                                                                                                                                                            |                                                                                                                                                                                                                                                                                 |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The Altrix system is suitable for use in the electromagnetic environment specified below. The customer or the user of Altrix should assure that it is used in such an environment. |                                                                                                                                                                            |                                                                                                                                                                            |                                                                                                                                                                                                                                                                                 |
| Immunity test                                                                                                                                                                      | IEC 60601 test level                                                                                                                                                       | Compliance level                                                                                                                                                           | Electromagnetic environment-guidance                                                                                                                                                                                                                                            |
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                                                                                                     | $\pm 8$ kV contact<br>$\pm 15$ kV air                                                                                                                                      | $\pm 8$ kV contact<br>$\pm 15$ kV air                                                                                                                                      | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                                   |
| Electrostatic fast Transient/burst<br>IEC 61000-4-4                                                                                                                                | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for input/output lines                                                                                                     | $\pm 2$ kV for power supply lines<br>$\pm 1$ kV for input/output lines                                                                                                     | Main power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                              |
| Surge<br>IEC 61000-4-5                                                                                                                                                             | $\pm 1$ kV line(s) to line(s)<br>$\pm 2$ kV line(s) to earth                                                                                                               | $\pm 1$ kV line(s) to line(s)<br>$\pm 2$ kV line(s) to earth                                                                                                               | Main power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                              |
| Voltage dips, voltage variations and short interruptions on power supply input lines<br>IEC 61000-4-11                                                                             | 0% $U_T$ for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°<br>0% $U_T$ for 1 cycle<br>70% $U_T$ (30% dip in $U_T$ ) for 25 cycles<br>0% $U_T$ for 250 cycles | 0% $U_T$ for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°<br>0% $U_T$ for 1 cycle<br>70% $U_T$ (30% dip in $U_T$ ) for 25 cycles<br>0% $U_T$ for 250 cycles | Main power quality should be that of a typical commercial or hospital environment. If the user of the Altrix system requires continued operation during power main interruptions, it is recommended that the device be powered from an uninterrupted power supply or a battery. |
| Power frequency (50/60Hz) magnetic field<br>IEC 61000-4-8                                                                                                                          | 30 A/m                                                                                                                                                                     | 30 A/m                                                                                                                                                                     | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                                       |
| <b>Note:</b> $U_T$ is the a.c. mains voltage before applications of the test level.                                                                                                |                                                                                                                                                                            |                                                                                                                                                                            |                                                                                                                                                                                                                                                                                 |

## Guidance and manufacturer's declaration - electromagnetic immunity

|                                                                            |                                                                                                 |                                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
|----------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|-----------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <p>Conducted RF<br/>IEC 61000-4-6</p> <p>Radiated RF<br/>IEC 61000-4-3</p> | <p>3 Vrms<br/>6 Vrms in ISM bands<br/>150 kHz to 80 MHz</p> <p>10 V/m<br/>80 MHz to 2.7 GHz</p> | <p>3 V<br/>6 V in ISM bands</p> <p>10 V/m</p> | <p>Portable and mobile RF communications equipment should be used no closer to any part of the Altrix system, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter.</p> <p>Recommended separation distance</p> <p><math>D=(0.35) (\sqrt{P})</math><br/>80 MHz to 800 MHz</p> <p><math>D=(0.70) (\sqrt{P})</math><br/>800 MHz to 2.7 GHz</p> <p>where <math>P</math> is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and <math>d</math> is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey,<sup>a</sup> should be less than the compliance level in each frequency range.<sup>b</sup></p> <p>Interference may occur in the vicinity of equipment marked with the following symbol:</p>  |
|----------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|-----------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

**Note - 1-** At 80 MHz and 800 MHz, the higher frequency range applies.

**Note - 2-** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**Note - 3-** The ISM (industrial, scientific, and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz; and 40,66 MHz to 40,70 MHz.

<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast, and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Altrix system is used exceeds the applicable RF compliance level above, the Altrix system should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the Altrix system.

<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths are less than 3 V/m.



3T Medical or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: **Altrix, Clik-Tite, Mul-T-Blanket, Rapr-Round, 3T Medical**. All other trademarks are trademarks of their respective owners or holders.



3T Medical Systems, Inc.  
236 Fernwood Avenue  
Edison, NJ 08837  
USA